
CHROMagarTM Pasteurella

Click below:

• CHROMagar™
Pasteurella

For detection of Pasteurellaceae



Plate Reading

- Pasteurellaceae
→ Pink to mauve

- *E. coli*, coliforms
→ inhibited or blue to steel blue

For detection of Pasteurellaceae

Background

In the veterinary field, respiratory infections in cattle herds can cause stunted growth, decline in milk production and death of animals, leading to significant economic losses.

Bacteria of the Pasteurellaceae family are commensal to the upper respiratory tract of many vertebrate species. Among them, *Histophilus somni*, *Pasteurella multocida* and *Mannheimia haemolytica* are among the main bacteria associated with the pathogenic complex of bovine respiratory diseases. During infection, these species cause complications which can lead to sepsis and death of the animal.

CHROMagar™ Pasteurella was therefore developed to improve the detection of Pasteurellaceae from a bovine respiratory sample (eg. nasal swab, trans-tracheal aspiration (TTA) and lung). This chromogenic medium to aid in qualitative diagnosis allows the detection and isolation of Pasteurellaceae colonies by inhibition or differentiation of the annex flora.

CHROMagar™ Pasteurella medium can also be used for samples from the respiratory tract of mammals such as cats and dogs.

Medium Performance

1 UNIQUE CHROMOGENIC MEDIA

First and only available specific medium for Pasteurellaceae in the market.

2 LOWER WORKLOAD

Contrary to general media like blood agar, CHROMagar™ Pasteurella differentiates the target group from the other bacteria by the a distinctive mauve colony color.

3 HIGHLY SELECTIVE

Most of the environmental flora will be inhibited.

4 EARLY PLATE READING

Most isolates will have a noticeable colony size after just overnight incubation.

Medium Description

Powder Base	Total 41.2 g/L Agar 15.0 Salt 5.0 Peptones 21.0 Chromogenic and selective mix 0.2 Storage at 15/30 °C - pH: 7.0 ± 0.2. Shelf Life > 18 months
Supplement S1 (Included in the pack)	Liquid form 8 mL/L Storage at 15/30 °C Shelf Life > 18 months
Supplement S2 (Included in the pack)	Powder form 1 g/L Storage at 2/8 °C Shelf Life > 18 months

Usual Samples	nasal swab, trans-tracheal aspiration (TTA), bronchial lavages, lung.
Procedure	Isolation & spreading method. Incubation 18- 24 h, 37 °C. Aerobic conditions.

Order References

Please use these product references when contacting your local distributor:
5000 mL pack PR012 (included in this reference: powder base PR012(B)
+ suppl. PR012(S1) + suppl. PR012(S2))

Manufacturer: CHROMagar

4 place du 18 juin 1940 75006 Paris - France
Email: CHROMagar@CHROMagar.com
Website: www.CHROMagar.com

Find your nearest distributor on
www.CHROMagar.com/contact

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