# Optitrol NAT SARS-CoV-2

en

## INTENDED USE

Optitrol NAT SARS-CoV-2 is a quality control (QC) sample intended to monitor the performance of in-vitro assays for determination of coronavirus of type SARS-CoV-2 RNA in human matrices. This QC was manufactured with the intention to always produce a positive result in NAT assays for SARS-CoV-2, which is helpful in determining the precision of testing systems and in identifying sources of variation. Optitrol NAT SARS-CoV-2 is an unassayed extraction and amplification control without target values.

Examples of possible reactivity (no target values) for this control can be found at <http://www.diamex.com>. Users of Optitrol QC samples have access to EDCNet, an Internet-based QC results monitoring system at <http://www.nrlquality.org.au/qconnect>.

The Optitrol NAT SARS-CoV-2 control is intended for use with assays from following companies:

Roche, R-Biopharm, altona Diagnostics, Qiagen, CerTest, Sansure Biotech, MyBioSource, CDC and others. Updates you will find at [www.diamex.com](http://www.diamex.com).

## CONTENTS

Optitrol NAT SARS-CoV-2 contains cell culture supernatants with intact SARS-CoV-2 viruses that were diluted in Viral Transport Medium (VTM). The cell culture supernatants with SARS-CoV-2 viruses were inactivated.

## STORAGE AND STABILITY

Optitrol NAT SARS-CoV-2:

- Can be used until the expiration date when stored at  $\leq -20\text{ }^{\circ}\text{C}$  (12 months from date of manufacture);
- Is stable for 5 days once opened providing it is closed tightly after use and returned to storage at  $2\text{--}8\text{ }^{\circ}\text{C}$ ;
- Do not dilute;
- Should always be stored upright;
- Must not be used beyond the expiration date.

## INSTRUCTIONS FOR USE

- Thaw at room temperature and use immediately or store at  $2\text{ }^{\circ}\text{C}$  to  $8\text{ }^{\circ}\text{C}$ ;
- Mix contents with shortly by vortexing before opening;
- After use, store at  $\leq -20\text{ }^{\circ}\text{C}$  or at  $2\text{ }^{\circ}\text{C}$  to  $8\text{ }^{\circ}\text{C}$ ;
- Optitrol NAT SARS-CoV-2 shall be treated as a patient sample, including the extraction and amplification steps;
- When using this control, the instructions of the test manufacturer must be followed exactly;
- The expected values when using the control must be determined by the user for the respective assay.

## LIMITATIONS

- Do not use this product where there is evidence of microbial contamination and/or high turbidity;
- Deviation from the recommended procedure may lead to unreliable results;
- This product is provided for quality assurance only and should not be used for calibration purposes;
- Possible causes of variation in Optitrol QC performance:
  - Different assay reagent lots;
  - Different calibrator lots used;
  - Different analysers used;
  - Different operators;
  - Time of day when testing.

## PRECAUTIONS AND DISPOSAL INFORMATION

- Optitrol NAT SARS-CoV-2 is for in-vitro diagnostic use only and should be used by trained personnel only;
- Optitrol NAT SARS-CoV-2 contains materials of human origin and should be considered potentially infectious. Observe universal precautions and handle this control like a patient sample;
- Wear protective gloves/ protective clothing;
- Do not pipette by mouth;
- Discard all materials in a safe manner that is in compliance with local and national regulations and in accordance to laboratory procedure for discard of clinical waste.

## WARNING

H226 Flammable liquid and vapour.

H302+H312+H332 Harmful if swallowed, in contact with skin or if inhaled.

## Prevention:

P210 Keep away from heat/sparks/open flames/hot surfaces. — No smoking.

P260 Do not breathe dust/fume/gas/mist/vapours/spray.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

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