


Zika Virus ViraStripe®

IgG und IgM Test Kits

The Zika Virus ViraStripe® IgG, IgM is an immunoassay employing different Zika Virus antigens. Relevant antigens are striped onto the membrane at defined positions. Specific reactivities at multiple banding positions yield

clear results to aid in the diagnosis of Zika Virus infections. Integrated negative control band, serum control, conjugate control and cut off control on each strip allow a safe, easy and reliable interpretation.



Easy interpretation through defined band positions

Economic - for low and high throughputs

Integrated cut off control on each strip

Native and recombinant antigens for excellent sensitivity and specificity

Envelope antigens: E, EIII, EIII*

Non-Structural Protein 1: important diagnostic marker

Fully automated strip processing available

Standardized evaluation with ViraScan® software

Zika Virus ViraStripe®:

The Standard for Specific Antibody Detection

Automated Evaluation with ViraCam® and ViraScan®



Order No.:
Zika Virus ViraStripe® IgG V-ZISGOK 50 tests
Zika Virus ViraStripe® IgM V-ZISMOK 50 tests

Order No.:
Positive control for Zika Virus ViraStripe® IgG V-ZISGPK 110µl
Positive control for Zika Virus ViraStripe® IgM V-ZISMPK 110µl
Negative control for Zika Virus ViraStripe® IgG/IgM V-ZISPKN 330µl