

Seronorm™ Immunoassay

The most complete routine immunoassay control in the market



Product features

- The "All-in-one" immunoassay control, combining hormones, cardiac and tumour markers
- Comprehensive product documentation, providing assigned values for more than 50 components from the most frequently used instruments on the market
- Available both as freeze-dried and liquid controls
- Three (freeze-dried)/four (liquid) clinically relevant levels – each level available separately
- Seronorm™ Immunoassay Liq low – a product which complements Seronorm™ Immunoassay with clinically relevant low levels of Cortisol, Ferritin, IgE, free T4, TSH, 25(OH) Vitamin D and Vitamin B12
- Long shelf life
- No preservatives or stabilizers added

Stability

- *Liquid material:*
 - 3 years stability when stored at $\leq -20\text{ °C}$ (with a few exceptions)
 - Opened vials (with a few exceptions):
1 month at $\leq -20\text{ °C}$ when refrozen within 30 minutes
10 days at $2-8\text{ °C}$
- *Freeze-dried material:*
 - 4 years stability when stored at $2-8\text{ °C}$
 - After reconstitution (with a few exceptions):
1 month at $\leq -20\text{ °C}$
10 days at $2-8\text{ °C}$

Components:

17- α -OH-Progesterone	IgE*
25(OH) Vitamin D*	IGF1
AFP	Insulin
Aldosterone	LH
AMH	Methylmalonic acid
Androstendione	Myoglobin
Anti-Tg	NT-proBNP
Anti-TPO	Ostase
β 2-microglobulin	Progesterone
β -hCG, total	Prolactin
CA 125	Prostatic acid phosphatase
CA 15-3	PSA, free
CA 19-9	PSA, total*
CEA	PTH, intact
CK-MB	SHBG
Cortisol*	T3, free
C-peptide	T3, total
DHEA-Sulfate	T4, free*
Digoxin	T4, total*
EPO	TBG
Estradiol	Testosterone
Ferritin*	Theophylline
Folate	Thyreoglobulin
FSH	Troponin I
hCG, total	Troponin T
hGH	TSH*
Homocysteine	Vitamin B12*

Refer to the package insert of currently available lots for specific analyte and stability claims.

* Also available in Seronorm™ Immunoassay Liq Low in clinically relevant low levels.

Order Information:

Art. no	Product Name	Size
207005	Seronorm™ Immunoassay Liq L-1	12 x 3 mL
207105	Seronorm™ Immunoassay Liq L-2	12 x 3 mL
207205	Seronorm™ Immunoassay Liq L-3	12 x 3 mL
207305	Seronorm™ Immunoassay Liq low	12 x 3 mL
206005	Seronorm™ Immunoassay Lyo L-1	12 x 3 mL
206105	Seronorm™ Immunoassay Lyo L-2	12 x 3 mL
206205	Seronorm™ Immunoassay Lyo L-3	12 x 3 mL