

With secured results against COVID-19

ViroQ SARS-CoV-2 – Real time PCR Kit

Your benefits

One-step multiplex real time SARS-CoV-2 Test

The ViroQ SARS-CoV-2 Kit is designed for specific detection of viral RNA Gene E and RdRP of SARS-CoV-2 in respiratory samples

Be confident of negative results

- False negative results due to a lack of cell material in the sample are excluded by simultaneous detection of a universally expressed human gene
- A mix of human RNA and SARS-CoV-2 reference RNA is used as positive control (included in the kit)

Our product features – Your benefits

- Highly specific and exclusive detection of the SARS-CoV-2 RdRP gene and Beta-CoV E gene
- One-Step real time PCR assay – easy to set up
- Results in 90 minutes – save time
- Works with common RNA-Isolation kits – be independent from single suppliers
- Ambient shipping condition – no extra costs and handling of dry ice
- Quality made in Germany – high supply chain security

Compatible with most real time cyclers


- The multiplex real time PCR kit uses three colours that are compatible with most common real time cyclers on the market
- Tested on QuantStudio 6 (Thermo Fischer Scientific), CFX 96™ (Bio-Rad) Lightcycler® 480 (Roche)

Workflow and Interpretation



Example for positive result

RdRP gene: FAM
E gene: CFO560 / HEX / VIC / JOE
Internal control: CFR610 / Texas Red / ROX



Ordering Information

Cat. No.	Kit	Size
728250	ViroQ SARS-CoV-2	96 Tests
728251	ViroQ SARS-CoV-2	480 Tests



BAG Diagnostics GmbH
Amtsgerichtsstr. 1-5
35423 Lich / Germany

Ordering:

Phone: +49 (0) 6404 925 450
Fax: +49 (0) 6404 925 460
Email: order@bag-diagnostics.com

Customer Service:

Phone: +49 (0) 6404 925 125
Fax: +49 (0) 6404 925 421
Email: info@bag-diagnostics.com